HUMAN CARE

Your Life. Your Way.

User manual / Bruksanvisning /
Notice d'utilisation / Benutzerhandbuch /
Handleiding

Slings / Selar / Harnais / Gurte / Tilbanden
Toileting Sling 25080 | 25085, Toileting Sling High 25090

Content

ENGLISH

User Manual	3
Safety Instructions	4
General Description / Intended Use	5
Technical Specifications	6
Product Label	7
Compliances and Standards	8
Suspension Alternatives	8
Hanger Bars for Human Care Lifts	9
Compatibility	10
Sling Application	11
Troubleshooting	14
Care Instructions	14

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User Manual

Important!

You must read the user manual for your device prior to usage. Keep this booklet and information handy for future use.

Symbols

WARNING symbol indicates a potentially hazardous situation which, if not avoided, could result in serious injury to the user or damage to property and/or the device itself.

INFORMATION symbol indicates recommendations and information for proper, trouble-free usage.

WEIGHT CAPACITY symbol indicates maximum user weight for the product. Do not exceed this weight in any circumstance.

CE MARKING symbol indicates product meets requirements of EU Directive 93/42/EEC (MDD) and/or EU Regulation 2017/745 (MDR).

MEDICAL DEVICE symbol indicates that the product is a medical device according to the definition of medical device in EU Directive 93/42/EEC (MDD) and/or EU Regulation 2017/745 (MDR).

RECYCLE according to local regulations.

READ INSTRUCTIONS and ensure you fully understand them before using this product.

Product modification is not permitted.
Before use, always check that the product is not damaged or worn.

Human Care's products are continuously being developed and updated to provide our customers with the highest quality. We reserve the right to make design changes without prior notice.

Always make sure that you have the most recent version of the manual which is available for downloading from our website at www.humancaregroup.com.

Contact your local distributor if you have any questions about the product and its use.

Notice to user/patient in case of serious incident

Any serious incident that has occurred in relation to the product, should be reported to the local contact, who reports to the manufacturer, and the competent authority of the country in which the user/patient is established.


Safety Instructions


The process of lifting and moving a person always involves risks. Therefore, you must read the user manuals for the lift and the lifting equipment carefully. Always make sure that the lifting equipment you are using is intended for use with the lift. As a care provider, you are responsible for the client's safety and you must know whether the client will be able to cope with the lifting process. The sling must be fitted to meet the client's specific needs. If you are in any doubt, please contact the supplier.

- Carry out a careful visual inspection of the sling every time you use it.
- Check that the right sling is being used for the client. The model, size and material must meet the client's needs.

- Plan the moving and lifting process so that it is as safe as possible.
- Before raising the client, check that the sling is correctly attached to the lift/hanger bar once the lift straps are stretched.
- It is important that the hooks of the hanger bar are in the correct position, with the opening of the hook facing upwards on both sides to securely hold the straps in place. Recheck the position of the hooks once the lift straps are stretched.
- It is important to always check that the loops are at the bottom of all hooks (see pictures) and securely placed below the locking pin on all hooks.

Correct

Incorrect

- Check that the client is sitting securely in the sling before starting to lift him/her.
- Check that the lift is balanced when moving, raising and/or lowering the client. If you are using a Human Care lift with telescopic hanger bar, it is important that the bar hooks are in the same position on both sides.
- Our slings have been designed and tested for indoor use with one client at a time and must only be used for transferring a person. This product is not a swing.
- Make sure no sharp objects are in contact with the sling.

General Description / Intended Use

Toileting sling

The Human Care Toileting Sling 25080 and the Human Care Toileting Sling Net 25085 are used primarily to take patients to the toilet. The slings have a reinforced back section and padding under the arms. The Toileting slings can also be used to move patients to wheel chairs and moulded seats.

The slings can be combined with the Human Care waist belt 25145/25146. The Toileting Sling 25080 and 25085 are designed for two-, three-point suspension.

Toileting Sling High

The Human Care Toileting Sling High 25090 is a sling with neck support specially designed for taking patients to the toilet. The sling has a reinforced back section and padding under the arms. It can also be used to move patients to narrow wheelchairs and moulded seats and in the bath and shower.

The 25090 model is made from low-friction polyester mesh. The advantage of this material is that the sling is easy to fit and dries quickly after washing or after the patient has taken a bath. The sling can be used with the Human Care waist belt 25145/25146. The Toileting Sling High 25090 is designed for two-, three- or four-point suspension.

Toileting Sling and Waist Belt

Both Toileting Slings can be completed with Human Care's Waist belt, with or without buckle, by threading the belt through the straps on the inside of the sling before application. When the user is sitting in the sling, close the belt (waist belt without buckle art. no.: 21545, alternatively with buckle art. no.: 25146).

A kit is available which includes a belt with buckle and sling. See below for kit specification.

Toileting Sling & Toileting Sling High Net w. belt


Product	Art. Nr.	Size	Consists of art. nrs. (waist belt and sling)	
Toileting Sling w. belt, JS	25080B-JS	Junior small	25146-40B	25080-JS
Toileting Sling w. belt, JM	25080B-JM	Junior medium	25146-60B	25080-JM
Toileting Sling w. belt, S	25080B-S	Small	25146-80	25080-S
Toileting Sling w. belt, M	25080B-M	Medium	25146-100	25080-M
Toileting Sling w. belt, L	25080B-L	Large	25146-110	25080-L
Toileting Sling w. belt, XL	25080B-XL	Extra large	25146-130	25080-XL
Toileting Sling High Net w. belt, JM	25090B-JM	Junior medium	25146-60B	25090-JM
Toileting Sling High Net w. belt, S	25090B-S	Small	25146-80	25090-S
Toileting Sling High Net w. belt, M	25090B-M	Medium	25146-100	25090-M
Toileting Sling High Net w. belt, L	25090B-L	Large	25146-110	25090-L
Toileting Sling High Net w. belt, XL	25090B-XL	Extra large	25146-130	25090-XL

Slings


Technical Specifications

It is crucial to use the right size of sling for the client. To help differentiate our sizes, they are colour coded for quick identification. The colour code is found next to the sling's label in the shape of a loop. This loop can also be used to properly store or hang the sling.

Toileting Sling

Toileting Sling High

Item	Part No.:	Size	Colour code	A (cm/ inch)	B (cm/ inch)	C (cm/ inch)	Weight Capacity* (kg/lbs)	Weight Capacity (kg/lbs)
Toileting Sling	25080-JS	Junior small	Purple	28/11	62/24,4	22/8,7	15-20/33-44	300/661
Toileting Sling	25080-JM	Junior medium	Grey	28/11	67/26,4	26/10,2	20-35/44-77	300/661
Toileting Sling	25080-S	Small	Red	32/12,6	74/29,1	30/11,8	30-60/66-132	300/661
Toileting Sling	25080-M	Medium	Yellow	48/18,9	90/35,4	33/13	55-90/121-198	300/661
Toileting Sling	25080-L	Large	Green	49/19,3	97/38,2	33/13	85-120/187-265	300/661
Toileting Sling	25080-XL	Extra large	Blue	56/22	100/39,4	38/15	110-140/243-309	300/661
Toileting Sling Net	25085-JM	Junior medium	Grey	28/11	67/26,4	26/10,2	20-35/44-77	300/661
Toileting Sling Net	25085-S	Small	Red	32/12,6	74/29,1	30/11,8	30-60/66-132	300/661
Toileting Sling Net	25085-M	Medium	Yellow	48/18,9	90/35,4	33/13	55-90/121-198	300/661
Toileting Sling Net	25085-L	Large	Green	49/19,3	97/38,2	33/13	85-120/187-265	300/661
Toileting Sling Net	25085-XL	Extra large	Blue	56/22	100/39,4	38/15	110-140/243-309	300/661
Toileting Sling High Net	25090-JM	Junior medium	Grey	28/11	73/28,7	70/27,6	20-35/44-77	300/661
Toileting Sling High Net	25090-S	Small	Red	32/12,6	82/32,3	74/29,1	30-60/66-132	300/661
Toileting Sling High Net	25090-M	Medium	Yellow	48/18,9	92/36,2	83/32,7	55-90/121-198	300/661
Toileting Sling High Net	25090-L	Large	Green	49/19,3	99/39	90/35,4	85-120/187-265	300/661
Toileting Sling High Net	25090-XL	Extra large	Blue	56/22	101/39,8	100/39,4	110-140/243-309	300/661

**The weights given are only guidelines and are not a substitute for individual tests. Discrepancies can occur, tolerances ± 1 cm, 0,3".*

Product Label

The label on the sling contains the following information:

- Human Care company logo
- This way up arrow and marking the centre line of the sling
- CE and MD marking
- Product name and size
- Picture showing the sling in use
- Part number
- Safe Working Load
- “This side facing out” – outside of the sling, not in touch with the client
- Manufacturer
- Date and country of production
- Instruction labels
 - Do not use fabric softener
 - Visual inspection
 - Latex Free
 - Read the manual before use
- Serial number
- Washing instructions
 - Machine Wash Temperature
 - Do Not Bleach
 - Tumble Dry, Low Heat
 - Do Not Dry Clean
 - Never Iron
- Owner name. (For example, write the client’s name or another unique identifier)
- Indicate bi-annual inspection
- UDI: (01) EAN/GS1/GTIN (11) manufacturing date (21) serial number
- Unique Product Code
- QR code to user manual

HUMAN CARE
The Life You Care For

UP

CE
MD

Toileting Sling High Net, Medium

REF 25090-M

SWL < 300kg / < 661lbs

This side facing out

Human Care HC AB
Arstaängsvägen 21B
117 43 Stockholm
Sweden
<http://www.humancaregroup.com>

2023-01-30

SN 335629220013

75°C / 167°F

#

Due date next periodic inspection (every 6 months), mark:

2023	1	2	3	4	5	6	7	8	9	10	11	12
2024	1	2	3	4	5	6	7	8	9	10	11	12
2025	1	2	3	4	5	6	7	8	9	10	11	12
2026	1	2	3	4	5	6	7	8	9	10	11	12

UDI (01) 0 7350011 05506 7 (11) 220130 (21) 335629220013

UPC 0816060801325

Compliances and Standards

Human Care is an ISO 13485:2016 certified Swedish medical device company. The Quality Management System is in compliance with US 21 CFR part 820.

The product is CE marked in accordance with EU

Directive 93/42/EEC (MDD) and/or EU Regulation 2017/745 (MDR), as class I medical device.

The CE mark is on the product.

This sling has been tested in accordance with ISO 10535.

Suspension Alternatives

Our slings are designed with the premise that all people are different with varying lifting requirements, needs or desires.

To make all our clients feel as comfortable as possible in our slings, most of them can be used with 2 or 3-point suspensions. Some slings can also be used with a 4-point suspension. Below the options are displayed.

2-point suspension

Suitable for most lifting situations providing a good seating position. 2-point suspension is available with floor lifts, fixed and portable ceiling lifts.

Two-point suspension with crossed leg support: Thread the strap on one leg section through the strap on the other leg section before attaching them to the hanger bar.

Two-point suspension with overlapping leg support:

This method is suitable for lifting clients with amputated legs, for example. Extension straps can be attached to the lifting straps to prevent the client from leaning too far back in the seated position.

Two-point suspension with separate leg support:

Put the client's legs in each leg section and attach the lifting straps to the hanger bar. Note: This method increases the risk of the client sliding out of the sling



Two-point suspension with Y-connection of the leg parts:

To create, for example, more space around the groin, change the pressure under the user's thighs or adjust the user's leg in relation to the hip, a Y-connection of the leg parts can be used.

Thread the left leg section through the leg support loop that is sewn on the outside of the right leg section, so that the leg parts go parallel to each other. Then attach the left leg part to the left hook of the hanger bar. Then do the same with the right leg part that is attached to the right hook of the hanger bar. This creates a Y connection (the connection looks like a Y from the front, see illustration).


An individual assessment of the client's needs must always be carried out before choosing an alternative connection method.

3-point suspension


Offers an elongated sitting position and more even weight distribution. The 3-point suspension is suitable especially for clients carrying an increased risk of pressure sores or who need additional space for the upper body.

Thread the left leg section through the loop on the right leg section. Attach the leg section lifting straps to the leg hooks and the back section lifting straps to the lift.

4-point suspension

A spacious sitting position making it suitable for larger clients or when using a stretcher. 4-point suspension is available with both fixed and portable ceiling lifts or floor lifts.


Suspension with belt


The sling can be completed with Human Care's Waist belt art.nr: 21545 or alternatively art nr: 25146 by threading the belt through the straps on the inside of the sling. When the user is sitting in the sling, close the belt round the users waist. When the belt with buckle is closed, a click sound indicates that the buckle is closed.


Hanger Bars for Human Care Lifts

HeliQ
(7275x)

Roomer S / Altair
(55751 45 cm, 17,7" / 55764 55 cm, 21,7")

HeliQ
(92938)

Accessory	Description	Part no.
2-point hanger bar	For stationary lift HeliQ. 3 widths:	35 cm/13,8": 72754 45 cm/17,7": 72755 60 cm/23,6": 72757
2-point hanger bar	For portable lifts Roomer S / Altair, 2 widths:	45 cm/17,7": 55751 55 cm/21,7": 55764
4-point hanger bar	For stationary lift HeliQ	92938
4-point suspension	4-point suspension accessory (49 cm). Enables 4-point lifting with the Altair and Roomer S lifts. The suspension bar is placed on the lift's single leg hook.	55619

Compatibility

Human Care Group slings are designed for use with all Human Care lifts and hanger bars. The products are also compatible with products of other suppliers utilizing the same attachment method for securing slings to a hanger bar and lift. However, due to the wide range of products worldwide, Human Care cannot be responsible for any errors that may result from improper use or application of combinations of slings, hanger bars and lifts of other suppliers. The use of the combined system is at the prescriber's own risk and liability.

To ensure the safety of the client and caregiver, Human Care recommends the following minimum requirements:

1. When using slings with lifts and hanger bars, made by manufacturers other than Human Care, the combined system needs to be designed with the same sling attachment method for the hanger bar of floor lifts and ceiling lifts, meaning loop style slings are combined with a loop style hanger bar, clip style slings are combined with a clip style hanger bar.
2. The involved patient lifts and sling must be CE marked in accordance with the EU Directive 93/42/EEC and/or EU Regulation 2017/745 and comply with standard

requirements for patient lifts as stated in EN ISO 10535.


3. The caregiver must follow the manufacturer's instruction manuals and recommendations with respect to the use, care, client, caregiver, cleaning, and inspection of the slings and lifts.
4. Different products on the assembled lift system: i.e., lift, hanger bar, sling, scale and other accessories may have different maximum loads. It is always the lowest maximum allowable load indicated for the respective products and the system that applies for the entire system. Always check the markings on all individual products in the system.
5. Individual risk assessment, including a practical test of the system and validation with the specific combination used are required to affirm that the size of the sling is correct for the intended use and compatible to the width and design of the hanger bar. The combined system needs to be documented by the competent assessor.

If there are any questions, please contact your local Human Care office or distributor.

Sling Application

Lifting the patient from a lying position


1. Lay the patient on his/her side. Position the sling so that the centre of the sling lines up with the patient's spine. Take note of the marking on the sling (the arrow pointing up, outside and the centre line). Roll the patient over onto his/her other side and pull the sling through.

Lay the patient on his/her back. If the bed is height-adjustable, it is a good idea to raise the head of the bed so that the patient is in a semi-sitting position to ease placing the sling behind the user's back.

2. These slings should preferably be used with two-point suspension. In this case, the loop on the right leg section is not used. Pull the leg sections through between the patient's legs. If the patient is being moved to the toilet, you can undo their trousers at this point. Cross the leg sections. Lower the lift to the height of the patient's stomach. Attach the lifting straps to the lift bar. Hold the lift away from the patient's face with one hand as you raise it above the patient's head.


Alternative: The Toileting Sling can also be used with three- or four-point suspension, but due to the big aperture it is only possible with patients who have good body posture and awareness for performing the transfer.


The arms of the patient are outside the sling during lifting to prevent sliding out of the sling.

Slings

3. Pause in this position and check that all the lifting straps are securely attached to the hooks.

4. Continue lifting the patient. It may be necessary to lower the bed slightly. Raise the patient's legs slightly while turning him/her away from the bed. If you are in the toilet, follow the instructions given in the section on lifting the patient from a sitting position.

Lifting the client from a sitting position

1. Pull the sling down behind the patient's back. Take note of the marking on the sling (the arrow pointing up, outside and the centre line). Pull the leg sections through to the front.


2. These slings should preferably be used with two-point suspension. In this case, the loop on the right leg section is not used. Pull the leg sections through between the patient's legs. If the patient is being moved to the toilet, you can undo their trousers at this point. Cross the leg sections. Lower the lift to the height of the patient's stomach. Attach the lifting straps to the lift bar. Hold the lift away from the patient's face with one hand as you raise it above the patient's head.

Alternative: The Toileting Sling can also be used with three- or four-point suspension, but due to the big aperture it is only possible with patients who have good body posture and awareness for performing the transfer.



3. Pause in this position and check that all the lifting straps are securely attached to the hooks.

4. Continue lifting the patient. It is a good idea to put your arm around the patient while you are moving him/her to create a greater sense of security.

If you are at the toilet, you can pull the patient's trousers down as far as they will go while the patient is still sitting in the sling. Then lower the patient onto the toilet, lower the lift to the height of the patient's stomach and disconnect the leg sections. (If the patient can sit without support, remove the sling completely.) Raise the lift above the patient's head. Pull the patient's trousers right down. The patient can now sit on the toilet supported by the sling.

The arms of the patient are outside the sling during lifting to prevent sliding out of the sling.

Troubleshooting

Description	Action
The client is at risk of sliding out	<ul style="list-style-type: none">• The sling may be too large.• The sling is not correctly fitted- the length of the loops should be stretched at the same time.
The client is unbalanced when lifted	<ul style="list-style-type: none">• The instructions on the individual test log have not been followed or a new assessment needs to be carried out.
The sling is too small	<ul style="list-style-type: none">• This can be uncomfortable for the client. The sling can cut into the client's flesh and/or his/her back may not have enough support.

Note: It is important that the right sling is used for each client and that it is individually tested to ensure that it meets the client's need and the transfer situation.

Care Instructions

Check for wear and tear

Check the sling for wear and tear. It is important to check the fabric, the edging, the seams and the lifting straps.

Damaged lifting equipment must not be used.

Cleaning

Human Care slings are intended for individual use. To avoid infections or transmission of bacteria, it is important to clean the slings properly. All slings should be machine washed and dried as soon as possible. Slings should always be cleaned after use in a wet situation.

Always refer to the label of the sling before washing.

Periodic inspection

A periodic inspection must be carried out at least every 6 months. See Human Care's periodic inspection log which can be found on our website at www.humancaregroup.com

Product lifetime

The product's expected lifetime 3 years or 150-180 times washing, if the product is used as intended and maintained according to the manufacturer's instructions, depending on the intensity of use and maximum load applied during use. If the product label is no longer legible, the product should be discarded.



Warranty and support

If you need information or support, please contact your Human Care office, www.humancaregroup.com, or distributor.

Recycling




Once the product has reached the end of its life, it should be disposed according to your country's rules, especially for the separate collection of electrical and electronic products.

Correct disposal helps prevent negative consequences for the environment and human health.

Overview of Sling Material Care		Polyester	Polyester Net	Trevira / Soft Air Mesh	Polyester/Polyester Net	Polyester/Cotton	Ventilated Polyester	Cotton	Polyester Tricot	Coated Nylon
	Machine Wash Maximum Temperature 60°C (140°F) (slings with padding)	●	●							●
	Machine Wash Maximum Temperature 75°C (167°F)	●	●	●	●		●		●	
	Tumble Dry Low, Maximum Temperature 30°C- 54°C (86°F- 130°F)									●
	Tumble Dry Low, Maximum Temperature 60°C (140°F)	●	●	●	●		●		●	
	Do Not Tumble Dry					●		●		
	Do Not Iron	●	●	●	●	●	●	●	●	●
	Do Not Dry Clean	●	●	●	●	●	●	●	●	●
	Do Not Bleach	●	●	●	●	●	●	●	●	●
	Do Not Use Fabric Softener	●	●	●	●	●	●	●	●	●
	Do Not Wash					●		●		
	Dispose of After Use					●		●		
	Disinfection Permitted (see list of approved disinfectants)	●	●		●		●		●	

Disinfection Procedure

Human Care slings are intended to be used with one patient at a time. To avoid possible transmission of infectious diseases, it is important to clean the sling if it has become soiled or if it is to be used by a different patient.

Approved Disinfectants	Solutions to Avoid
Chemical Disinfection: 70% Ethanol, 45% isopropanol or similar	 Steam sterilization
Thermal Disinfection: Wash 65°C (149°F) for a minimum of 10 minutes ¹	 Iodine based disinfectants
Thermal Disinfection: Wash 71°C (160°F) for a minimum of 3 minutes ¹	 Gas sterilization
Note: Add 8 minutes to both thermal disinfection options to ensure heat has penetrated the fabric. Finish with a rinse cycle. ¹	

¹ Ref: HSG (95) 18: Hospital Laundry Arrangements for Used and Infected Linen, Department of Health, NHS (National Health Service, England)

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